

Carestream

Focus 35/43C Detectors
with Image Suite V4 Software



Wireless. Digital. Affordable.

Upgrade Your Practice to DR Imaging.



SMART IMAGING SOLUTIONS

The Budget-Friendly Way to Upgrade to DR Imaging.

Stepping up to full digital X-ray is a tremendous way to speed your workflow, increase productivity and support a higher standard of care. But the cost of upgrading can be a real barrier.

That's why we've introduced our highly economical CARESTREAM Focus 35C/43C Detectors, powered by CARESTREAM Image Suite Software. It delivers the productivity power of wireless operation and the image quality only DR technology can provide.

Best of all, the Focus Detector slides right into your existing analog equipment – protecting your current investment.



Advantages at a Glance

- This is an affordable way to upgrade to full digital imaging.
- The superb image quality of DR supports greater diagnostic confidence.
- Wireless design boosts workflow.
- Image Suite V4 Software includes a full suite of veterinary measurement tools.
- Choose from two sizes: 35x43 and 43x43 cm.





Image Suite: Optimized for Veterinary Applications.



*Image Suite software adapts to the way you work –
with a custom-designed veterinary interface.*

Many digital-imaging systems marketed to vets come with a generic, human-based clinical interface that requires you to adapt. By contrast, Image Suite software adapts to the way you work – with a custom-designed veterinary interface.

Our order-entry screens have the fields you need, for selection of breed, species and more. You can also input unique identifiers like tattoos and chip IDs. Veterinary-specific reference images are provided to support accurate positioning, while sample veterinary X-rays provide guides for enhanced image quality. Our customized animal-image processing database fine-tunes processing parameters for specific species, breeds, body parts and views. Specialized veterinary measurement tools include:

- Hip Dysplasia Measurement Tool
- VD Thorax Measurement Tool
- Vertebral Scale Measurement Tool
- Tibial Tuberoses Advancement Measurement Tool
- Tibial Plateau Leveling Osteotomy Measurement Tool
- Laminitis Measurement Tool

The Power Behind the Detector: Image Suite Software.

The Image Suite software platform combines advanced image processing and robust functionality for maximum efficiency. Image Suite Software:

- Provides a simple and intuitive user interface custom-designed for veterinary practices.
- Allows automated workflow with uniform protocols.
- Adds extra Cybersecurity with the Windows 10 Platform.
- Offers an optional Mini-PACS module with a broad tool set, allowing superb image management and analysis.
- Performs backups to CDs and USB drives for reliable off-line storage.
- Offers a tablet-viewing option* for image review with patients to help you better explain your diagnosis and treatment planning.

* Tablet-Viewer Software is a purchasable option and requires a separate license.

Image Suite Software's Advanced Reading and Measurement Tools:

Hip Dysplasia



Tibial Tuberosity Advancement



Vertebral Heart Score



NO MATTER WHICH ROAD YOU'RE ON,
WE'LL GET YOU WHERE
YOU NEED TO GO.

SMART IMAGING SOLUTIONS



Say Goodbye to the Drawbacks of Analog Imaging.

Our Focus Detector with Image Suite V4 Software is the ideal way to tap the power of digital – and it's budget-friendly.

Additional Focus Detector advantages include:

- 35x43 cm and 43x43 cm sizes with Cesium Iodide scintillator for low-dose applications.
- Wireless design to help minimize the infection risks and trip hazards of cables.
- Beam Detect Mode to eliminate cable connection to the generator, alleviating OEM service concerns; image capture begins automatically when the detector senses the X-ray exposure.

Focus 35C/43C Detector Specifications:

Model Name	Focus 35C	Focus 43C
Size	35x43 cm	43x43 cm
Weight	3.43 kg/7.56 lbs.	4.68 kg/10.31 lbs.
Battery Life	6.5 hrs	6 hrs
Spatial Resolution	3.3 lp/mm	3.6 lp/mm
Scintillator	CsI	CsI



A Community of Service and Support.

For dependable service, look to our Customer Success Network. We work continuously to improve your imaging performance, help you to innovate as needs change, and make the most of your budget and resources. Carestream's Customer Success Network surrounds you with a dynamic team of experts, with a Single Point of Entry for easy, customized access to the right people in every situation. You and your patients will benefit from the expertise and best practices only Carestream can deliver.



carestream.com



Rx only

© Carestream Health, Inc., 2025. CARESTREAM is a trademark of Carestream Health. CAT 2000 338 05/25
Not all products or features are available in all countries. Please speak to a Carestream representative for more information.

Carestream