

Carestream

DRX-Rise Mobile X-ray System

# Digital. Powerful. Efficient.

Mobility Meets Versatility.



SMART IMAGING SOLUTIONS

# Mobile Imaging On The Rise.

## Elevate Your Imaging with Greater Mobility.

Looking for a mobile X-ray system to replace your existing unit or add to your current fleet? Is your facility still relying on a CR-equipped mobile unit?

The ideal solution would be the addition of a new, fully digital mobile X-ray system. And now it's easier and more affordable than ever to obtain.

The CARESTREAM DRX-Rise Mobile X-ray System is a fully integrated DR mobile X-ray solution that provides an easy workflow, elevates productivity and produces outstanding images.



### Advantages at a Glance

- Dual displays streamline workflow and boost productivity
- Quickly reposition the cart from the tube head or remote
- Provides superb maneuverability, even in tight and cluttered environments
- Offers seamless compatibility with a range of Carestream DR Detectors
- Utilizes user-friendly ImageView software

## Benefits Throughout Your Facility.

The DRX-Rise offers compelling features for every stakeholder in your imaging operation.



#### Rad Admin.

- A variety of detector options to align with their facilities' needs.
- Future-proof design to prevent technology obsolescence.



#### Radiographer

- Easy-to-manuever system allows for effortless navigation.
- Secondary drive controls allow quick adjustments while positioning.
- Dual displays and extended battery life boost productivity.



#### Radiologist

- Pristine images with a variety of image-processing options for unrivaled diagnostic confidence.
- A wide range of specialized viewing options for clear and precise imaging, whether in the ER, OR or ICU.



#### Patient

- High-quality X-ray images support accurate diagnosis and commencement of treatment.
- The convenience, comfort and minimal disruption of bedside imaging.

# Fully Integrated for Value and Performance.


The DRX-Rise was designed with input from radiographers to deliver improved workflow and seamless operations. And, its full and robust set of features offers a remarkable combination of productivity, power and convenience.

## Controls on the tube head

and the remote put system positioning right at your fingertips, allowing you to quickly reposition the cart.

## Dual Displays

Tube head display provides secondary workzone. Large main monitor makes it easy to navigate, review, and manipulate images.



**Remote exposure** via the handheld wireless remote exposure switch allows the radiographer to acquire exposures from anywhere in the area.

## In-bin detector charging

saves time and means the detector is always ready to go.

**Onboard storage** keeps all the supplies radiographers need close at hand.

## Functional LED lighting

indicates system status at a glance.

**ImageView software** utilizes our powerful and intuitive software platform to elevate workflow and deliver unparalleled image quality.

**Lithium battery** provides faster recharging and long life – up to 50 km of driving or 1,000 exposures per single charge.

**Compact design** and agile operation let you navigate throughout your facility with ease.



## Experience Powerful Software and Detector Compatibility.



### Harness the Power of ImageView Software.

The DRX-Rise is driven by our industry-leading ImageView Software. Its Eclipse engine applies proprietary algorithms and artificial intelligence to drive superb image quality and provide diagnostic confidence. Additional image-processing options include:

- Tube and Line Visualization
- Pediatric Image Optimization
- Pneumothorax Visualization
- Bone Suppression
- SmartGrid

### High Performance Runs in the Family.

When evaluating detector options, our family of detectors has the answer for your facility, whatever your patient population, workflow, dose requirements, exam types, or budget. And the DRX-Rise is compatible with every one of them.\*

**DRX Plus Detectors** (the original shareable detectors) feature the X-Factor – the ability to share the detector with any other DRX System.

**DRX Plus 2530C Detector**, with a 98-micron pixel pitch, delivers all the advantages of the DRX Plus in a smaller-format, pediatric-friendly design.

**Focus 35C and 43C Detectors** deliver a highly affordable, wireless DR solution while providing excellent image quality and an IP44 rating. **Focus HD 35 and HD 43 Detectors** have a 100-micron pixel pitch to provide higher resolution for improved image details.

\*All detectors may not be available in all regions.





NO MATTER WHICH ROAD YOU'RE ON,  
WE'LL GET YOU WHERE  
YOU NEED TO GO.

SMART IMAGING SOLUTIONS

### The Fastest Route to a Successful Imaging Future.

As you travel to the future of imaging, what's your next step? Perhaps you need to accelerate your X-ray workflow. Or add a new imaging solution to your facility to provide better patient care.

Wherever you're headed, we have the smart solutions you need: Full-digital imaging rooms and leading-edge mobile units. Wireless, shareable detectors and CR-to-DR retrofit kits. Affordable CR systems. And advanced digital imagers.

### Start Mapping Your Route Today.

From small clinics, urgent care facilities and specialty practices to large hospitals and imaging centers, every facility will find a Carestream solution to meet the most specific needs. Visit [carestream.com](http://carestream.com) and start your journey to greater productivity and a higher standard of patient care.

### More Solutions from Carestream.

Carestream also offers Non-Destructive Testing solutions to help ensure quality and safety in the aerospace, petrochemical and other industries, as well as Contract Manufacturing to bring innovative coated products to market.

### A Community of Service and Support.

For dependable service, look to our Customer Success Network. We work continuously to improve your imaging performance, help you to innovate as needs change, and make the most of your budget and resources. Carestream's Customer Success Network surrounds you with a dynamic team of experts, with a Single Point of Entry for easy, customized access to the right people in every situation. You and your patients will benefit from the expertise and best practices only Carestream can deliver.



### Eclipse: Harness the Power of Decades of Innovation.

Carestream's Eclipse is the engine behind our innovative imaging software. It uses AI technology and proprietary algorithms to significantly increase the value of the entire imaging chain, from capture to diagnosis.

[carestream.com/drx-rise](http://carestream.com/drx-rise)



"Rx only"

© Carestream Health, Inc., 2023. CARESTREAM is a trademark of Carestream Health. CAT 2000 299 02/23

Carestream