

Carestream

Vita Flex CR System

# Speed. Productivity. Quality.

Upgrade Your Practice from Film to Digital.



SMART IMAGING SOLUTIONS

# Big Performance for Smaller Facilities.

If yours is like many small practices or clinics, you realize that in-house digital imaging could provide you with tremendous value.

But the high cost of entry for digital may be holding you back – not to mention the limited floor space you have for installation of such a system.

Meet the CARESTREAM Vita Flex CR System – designed specifically for facilities like yours. It delivers high-quality digital images to speed workflow and improve patient care, and offers a range of throughput rates to match your workflow needs. Compact design makes it a great fit for smaller spaces. And its affordable price makes it possible on nearly any budget.

## Advantages at a Glance

- Eliminates the need for ongoing purchases of film and the use of harsh chemicals
- Saves hours compared to time-consuming film-based radiography
- Allows for the unit to be set up vertically or horizontally, with a unique space-saving design
- Works great for smaller facilities, because it's compact, lightweight and affordable





### FlexFit Means There Are Two Ways to Save Space.

Vita Flex's benefits for small offices don't end with its compact footprint. You can find the fit that's best for you. Vita Flex fits on a tabletop, desktop, on the floor or in the back of a van. You can position it in a traditional horizontal configuration, or set it up vertically, if that's a better fit for your space.

Now that's flexibility.



## Custom-Designed with Your Needs in Mind.

Why pay for capabilities you don't need and won't use? Large, complicated digital-imaging systems are overkill for your facility's needs. Carestream's Vita Flex CR delivers precisely – and only – the capabilities you require. Just consider these advantages:

### AFFORDABLE PURCHASE PRICE

- Go digital without breaking the bank.

### STREAMLINED OPERATIONS

- Say goodbye to ongoing consumable purchases.
- Forget the need to handle harsh chemicals.
- Now there's no more film to process.

### SMALL AND LIGHTWEIGHT

- Save valuable space with Vita Flex's compact footprint.
- Configure it vertically or horizontally to fit the tightest spaces.
- The system transports easily for mobile imaging.

### HIGH-QUALITY IMAGES

- Upgrade to superb digital images.
- Enjoy easy-to-use Image Suite acquisition software.
- Take advantage of powerful analysis tools in the software to improve diagnostics.

### SPEEDY PERFORMANCE

- View high-quality images in just minutes.
- Choose the throughput rate that best suits your workflow and exam types (30, 45 or 60 PPH).

### EXPERIENCE THE VERSATILITY

- Select from multiple cassette sizes.
- Handle all major exam types with confidence.

### OPTIONS

- The Dental Panoramic Kit allows fast and accurate dental imaging.
- The Long-length Imaging System enables flexible, high-quality long-bone image capture.

### MULTIPLE CR SCREENS AND CASSETTES

Ideal for use with the Vita Flex CR Systems, our line of CR cassettes and phosphor screens offers a wide selection and high performance to meet the needs of various CR exams and clinical applications, including long-length imaging, orthopaedics, pediatrics and more.

- 8 x 10 in.
- 10 x 12 in.
- 14 x 33 in.
- 14 x 14 in.
- 14 x 17 in.
- 24 x 30 cm
- 15 x 30 cm



# Flexibe and Reliable.

## Our Easy-to-Use Long-Length Imaging System.

Enjoy the benefits of digital capture and processing of long-length exams with our easy-to-use Long-Length Imaging System.

- Automatic image stitching
- Positioning flexibility
- Portability
- High image quality

### AUTOMATIC IMAGE STITCHING

Image processing algorithms are applied to optimize the image and the result is a single image. No manual alignment is required. The stitched images can be distributed to softcopy workstations or hardcopy laser imagers, just like any other scanned image.

### PORTABILITY

The wheeled caddy (optional) makes it easy to move the cassette to different locations and allows for imaging of supine long-length exams or upright exams.



### POSITIONING FLEXIBILITY

The vertical cassette holder (optional) allows capture of upright long-length exams. Exams can be captured as closely to one-half inch from the floor, which eliminates the need for patients to climb a step. The holder can be adjusted in one-inch increments (up to 74 inches from the floor) to accommodate each patient's height and body part to be imaged.

# Save Time. Save Money.

## Innovative Options for Installation and Service.

### SELF-INSTALL OPTION

- Set up and start scanning immediately with easy self-install.
- Go from carton to first image in less than 90 minutes.
- Or choose full-service installation from a certified technician.

### SELF-SERVICE OPTION

- Fast self-repair maximizes your uptime and budget.
- Fix it yourself with three easy, field-replaceable units.
- Or call for full service from a certified technician.

### IDEAL FOR MOBILE APPLICATIONS

Need to take Vita Flex on the road? Its small size and light weight make it a breeze to transport and use in the field.



[carestream.com](http://carestream.com)



"Rx only"

© Carestream Health, Inc., 2024. CARESTREAM is a trademark of Carestream Health. CAT 200 0053 02/24

# Carestream